

Your OneTouch® Verio®Sync Meter. Your iPhone®. Together in a few easy steps.

Follow these steps to sync your meter with an Apple® iPhone®, iPad®, or iPod touch®.



1. Get the OneTouch® Reveal™ mobile app.

- > Go to the App Store™ on your iPhone.
- > Type “OneTouch Reveal” in the search bar, tap **Search**, and locate the app.
- > Tap **Free**, then tap **Install app**.
- > The app icon will appear.
- > Open the app and accept the user license by tapping **I accept**.



2. Turn on your OneTouch® Verio®Sync Meter.

- > Press and hold the meter button (found on the right side of your meter).



Meter button

3. Connect your meter and your iPhone.

- > Activate *Bluetooth*® on your meter and your iPhone.
 - On your meter: press the meter button twice. This symbol  will appear on-screen.
 - On your iPhone: tap **Settings** and set Bluetooth to **On**. (For older operating systems, tap **General**, then tap **Settings**, then set Bluetooth to **On**.)
- > When Bluetooth is activated, your meter will appear in the iPhone **Devices** list. (The name will begin with “OneTouch.”)
 - If Bluetooth is already on, tap **Bluetooth** and the **Devices** list will appear.
 - Tap on the meter name and enter the PIN located on the back of your meter.



- Tap **Pair** and wait for **Connected** to appear next to the meter name.

4. Your first sync.

- > Open the OneTouch® Reveal™ mobile app on your iPhone and turn your meter on.
 - The sync symbol  will appear and blink on the meter screen. Your iPhone will display a **Syncing Data** message.
 - The app will display a message about meter time difference. Tap **Yes** to match the time on your meter and your iPhone.

You're all set!

Now you're ready to test and sync with your meter. Whenever you open the app on your iPhone, your blood sugar results will automatically appear in the app.†

*Meter must be on and in range (10 feet).

†Make sure Bluetooth is activated on both your meter and your iPhone.

The OneTouch® Verio®Sync Meter uses only OneTouch® Verio® Test Strips.

Don't have a prescription?

Contact your healthcare professional and pick up your OneTouch® Verio® Test Strips at the pharmacy.



The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by LifeScan, Inc. is under license. Other trademarks and trade names are those of their respective owners.

Apple, the Apple logo, iPad, iPhone, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. © 2013 LifeScan, Inc. Milpitas, CA 95035 6/13 DCF-13-2514A

Before using this product to test your blood sugar, please read the OneTouch® Verio®Sync Owner's Booklet. Have questions or need support? Visit www.OneTouch.com or call (888) 567-3003.

ONETOUCH®