

ONETOUCH  
Reveal®



## OneTouch Reveal® web app

### Helps you get more meaning out of your results

- Tracks your blood glucose, food, insulin dosing and activity data over time
- Suggests personalized tips to help you stay on track and in range\*
- Reports can be printed, ready for discussion at your next appointment with your healthcare professional
- Can be used on PC and Mac computers at [OneTouchReveal.com](http://OneTouchReveal.com)

\* Be sure to talk to your healthcare professional about the actions that are right for you.

# How to setup your OneTouch Reveal® web app

Before you begin, confirm that cookies are enabled in your browser.



Bluetooth

## Wireless Option

OR

## Cable Option

### 1. Connect wireless devices and activate your account

- Install the OneTouch Reveal® mobile app on your mobile device
- Open the mobile app and follow the on-screen instructions to setup your account and pair your OneTouch Verio Flex® meter
- Check your email inbox for a link to confirm your email address

For more details, refer to the OneTouch Reveal® mobile app **User Guide**.

### 2. Sign in to the OneTouch Reveal® web app on your computer

- Go to **OneTouchReveal.com** and enter the same **Username** and **Password** that you used in the OneTouch Reveal® mobile app
- Click **Sign In** and follow the on-screen instructions

### 1. Register and activate your account

- Go to **OneTouchReveal.com**, click on **Create a New Account** and follow the prompts
- Check your email inbox for a link to confirm your email address

### 2. Install Data Transfer Tool

- Unplug all OneTouch® USB cables from your computer
- **Sign In** to your OneTouch Reveal® account and click **help** from the top right of any screen
- Click the **Download** button under **Data Transfer Tool** and follow the prompts
- Select the appropriate devices
- Click **Finish**. After installation, you may be prompted to re-start your computer.

### 3. Transfer data from your OneTouch® meter

- Connect your meter to the computer with your OneTouch® USB cable (*cable sold separately*).
- Open the OneTouch Reveal® Data Transfer Tool and select **Detect Now**. When the device is found, a confirmation screen will appear.
- **Sign In** to the OneTouch Reveal® web application
- Click **Confirm** to transfer data. When the transfer completes, click **OK**



**You're all set! You can now access your data with the OneTouch Reveal® web app at OneTouchReveal.com**



## View your blood glucose readings

- Click the **Logbook** tab to see your blood glucose readings
- Click the **Patterns** tab to view a 14-day summary

**OneTouch Reveal**

**Patient Summary: 5/5/2016 - 5/18/2016**

Name: Patient 3, Kate    Patient ID: TTCP3    Overall Target: 70 - 180 mg/dL  
 DOB: 5/1/1987    Clinic: Tolland Training Clinic    Before Meal Target: 70 - 130 mg/dL  
 Gender: Female    Clinician: Clinic Man... TTU    After Meal Target: 70 - 180 mg/dL

Warning: Reports are intended for use by healthcare professionals with experience treating patients with diabetes only.

Use Patient Overall: 18 glucose readings with latest target (13.2%).  
 Use Patient Before Meal: 3 after meal glucose readings with latest target (before meal readings of 70).  
 Use Patient After Lunch: 100 glucose readings with latest target (28.2%).

Comprehensive Statistics

	Current Date Range	Previous Date Range	Change %
Blood Sugar (Mean)	104.0 mg/dL	103.3 mg/dL	0.6%
Standard Deviation	64.4 mg/dL	64.3 mg/dL	0.0%
Std. of Range	102	102	0%
High No. of Hypoglycemia day	1	1	0%
No. of Hypoglycemia Night	0	0	0%
After Meal Day	0	0	0%

Glucose by Time of Day

Logbook: 5/17/2016 - 5/30/2016

Name: Patient 3, Kate    Patient ID: TTCP3    Overall Target: 70 - 180 mg/dL  
 DOB: 5/1/1987    Clinic: Tolland Training Clinic    Before Meal Target: 70 - 130 mg/dL  
 Gender: Female    Clinician: Clinic Man... TTU    After Meal Target: 70 - 180 mg/dL

Warning: Reports are intended for use by healthcare professionals with experience treating patients with diabetes only.

Date	Overnight	Before Breakfast	After Breakfast	Before Lunch	After Lunch	Before Dinner	After Dinner	Bedtime
5/17/2016	94 @ 12:30 am 24 @ 1:31 am	113 @ 7:03 am		113 @ 11:30 am	210 @ 2:00 pm 99 @ 3:00 pm	83 @ 7:00 pm	73 @ 8:10 pm 96 @ 9:47 pm	102 @ 10:20 pm
5/18/2016		81 @ 7:00 am		114 @ 1:30 pm		105 @ 5:25 pm	96 @ 9:47 pm	138 @ 10:00 pm
5/19/2016	204 @ 14:58 am 140 @ 4:42 am	100 @ 6:40 am		112 @ 7:24 am		100 @ 5:25 pm	88 @ 9:30 pm 88 @ 10:49 pm	176 @ 10:58 pm
5/20/2016	51 @ 11:19 am	176 @ 6:27 am		87 @ 12:02 pm		146 @ 7:30 pm	144 @ 8:20 pm	176 @ 10:58 pm
5/21/2016	52 @ 12:47 am	65 @ 6:00 am		88 @ 12:53 pm		97 @ 5:53 pm	95 @ 7:24 pm	110 @ 9:25 pm
5/22/2016	62 @ 12:00 am	115 @ 6:30 am		63 @ 7:00 am		124 @ 4:30 pm	102 @ 7:12 pm	113 @ 9:25 pm
5/23/2016	52 @ 12:00 am	65 @ 6:00 am		62 @ 12:00 am		105 @ 5:00 pm	104 @ 8:20 pm	113 @ 9:25 pm
5/24/2016	47 @ 12:14 am	68 @ 7:17 am		155 @ 1:14 pm		179 @ 6:57 pm	144 @ 9:11 pm	204 @ 11:23 pm
5/25/2016	62 @ 12:00 am	65 @ 6:30 am		104 @ 6:51 am		104 @ 10:14 am	104 @ 11:23 pm	204 @ 11:23 pm
5/26/2016	62 @ 12:00 am	65 @ 6:30 am		100 @ 12:33 pm		78 @ 6:20 pm	78 @ 6:20 pm	102 @ 11:30 pm
5/27/2016	52 @ 12:00 am	65 @ 6:30 am		125 @ 12:00 pm		67 @ 6:10 am	67 @ 6:10 am	143 @ 10:31 pm
5/28/2016	52 @ 12:00 am	65 @ 6:30 am		103 @ 7:30 am		204 @ 6:00 pm	204 @ 6:00 pm	103 @ 11:19 pm
5/29/2016	52 @ 12:00 am	65 @ 6:30 am		102 @ 7:00 am		132 @ 6:12 pm	132 @ 6:12 pm	150 @ 10:35 pm
5/30/2016	43 @ 12:00 am	65 @ 6:45 am		111 @ 11:53 am		110 @ 6:45 am	306 @ 7:00 pm	140 @ 10:40 pm
5/31/2016	52 @ 12:00 am	65 @ 6:30 am		122 @ 12:00 pm		154 @ 6:14 pm	236 @ 7:00 pm	141 @ 11:14 pm
6/1/2016	52 @ 12:00 am	134 @ 10:30 am		177 @ 1:42 pm		83 @ 4:35 pm	114 @ 7:00 pm	114 @ 7:00 pm
6/2/2016	52 @ 12:00 am	65 @ 6:45 am		127 @ 8:48 am		112 @ 6:00 am	112 @ 6:00 am	101
6/3/2016	52 @ 12:00 am	65 @ 6:30 am		100 @ 12:40 pm		70 @ 5:30 pm	70 @ 5:30 pm	101



## Prepare for your next visit with your healthcare professional

- Click **Visiting your doctor?**  
The **Be prepared for your visit** screen will appear.
- Click **Download reports**  
To include reports, be sure to click the box next to **Include reports for your doctor** before you click **Download reports**.

If you'd like to change the list of reports you want to print for your visit, you can do so by clicking on **Settings**, then **Preferred Reports to Print** and then selecting the reports you need.

# Any questions? We're here to help



Visit **OneTouch.com**



Call us at 1 866 693-0599

Find the necessary system requirements at  
[OneTouchReveal.com/support/en\\_US/](https://OneTouchReveal.com/support/en_US/)

