

OneTouch Reveal® Mobile App

Compatible Devices and System Requirements

To use the OneTouch Reveal® mobile app you need:

- A compatible OneTouch® meter with *Bluetooth*® connectivity
- A compatible mobile device running a supported operating system
- To download the app to your mobile device from the App StoreSM or Google Play

Compatible OneTouch® Meters:

- OneTouch Verio Flex™ meter
- OneTouch Verio® Sync meter

Compatible Mobile Devices & Operating Systems:

- Apple iPhone®, iPad® or iPod touch®, running iOS 8 or later*
- Google Android™ devices with BTLE (Bluetooth® Low Energy) support and running Android Version 4.3, 4.4, 5.1 and 6.0.

In order to use the Data sharing (Share) feature within the OneTouch Reveal® mobile app, your mobile device will need Wi-Fi or cellular connectivity to send by SMS message or email. See your mobile provider tariff for charges on SMS or data usage.

*iPad2 does not have low-energy Bluetooth® and therefore will not pair with low-energy Bluetooth® meters or devices.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by LifeScan Scotland Ltd. is under license. Apple, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google, Google Play, Android and the Google Logo are registered trademarks of Google Inc.